



FEATURES

The **DESORB 5001**, packaged in convenient 5 Liters "pillows", is formed from polypropylene microfibers with a double-walled structure, which guarantees a long life and a remarkable resistance to tearing.

The product has a high resistance as it is unassailable by any chemical, toxic, acid, solvent etc.

The absorption capacity is very high and is therefore suitable to collect any harmful and dangerous chemical product.

Each pillow absorbs up to 4 liters of liquid without dripping and without losing the absorbed liquid when subjected to pressure.

CHEMICAL/PHYSICAL DATA

Composition:	polypropylene
Flammabaility:	440 °C
Absorbent power:	3.8-4 Lit / each.
Melting point:	150-200 °C
Volume:	5 lt
Weight:	550 gr
Dimensions:	40 x 40 cm
Packaging:	16 pz / box

USE

- Absorption of any chemical that escapes during tank loading and unloading operations.
- Absorption of any chemical that escapes from greasing or pickling plants.
- As a protective mat under annoying losses of chemical product, under cisterns, containers or baskets.

Insert the pillow under the source of purging or dispersions of valves, fittings and equipment that have losses.

DESORB 5001 prevents the chemical from dispersing on the work floor, keeping it clean and making the environment safer and more productive.

The absorbent starts working immediately and continues until it is completely saturated.

When the product is completely saturated, place it in a container and dispose of it according to the regulations for the pollutant.

